

AARON® 3250™

*Bovie*®

## ELECTROSURGERY

FOR THE HOSPITAL OPERATING ROOM AND SURGERY CENTER



ELECTROSURGICAL GENERATOR

### CUTTING EDGE

300 Watts of Cutting power

### CONVENIENCE

Ten (10) Doctor Presets for quick, uncomplicated setup between procedures

### VERSATILITY

Offers two (2) modes of Pure Cut, ten (10) modes of Blend, two (2) modes of Coagulation, Bipolar Coagulation and BovieFCFS™ technology allowing the Doctor to move almost effortlessly from Monopolar to Bipolar during procedure

### INTELLIGENCE

The Aaron 3250 is designed with BovieDFFS™ technology measuring tissue impedance 5,000 times a second for consistent energy delivery into varying tissue types

### SECURITY

The Aaron 3250 is designed with state-of-the-art safety features including Bovie NEM™ pad sensing technology which monitors the return electrode for optimum patient protection. BovieDED™ monitors all functions of the generator to ensure constant safety for the patient, doctor and OR staff





### RELIABILITY

The Aaron 3250 is backed by Bovie's two (2) year manufacturer's limited warranty

OTHER OPTIONAL  
ACCESSORIES

Disposable Pencils*	Disposable Grounding Pads*	Monopolar & Bipolar Footswitches
		
ESP1 & ESP7	ESRE & ESREC	BV-1253B & BV-1254B

Disposable Suction Coagulator*	Disposable Arthroscopic Electrode*	Resistick II™ Coated Blade*
		
SCH10 Handswitch	AR03	ES01T

Disposable and Reusable Electrodes*			
Blade	Ball	Needle	SuperCut™ Tungsten Needle
			
ES01	ES20	ES02	ES60R

\*More styles available

Industry-standard accessories fit most electrosurgical generators

Additional accessories available — contact your Bovie representative for details.

TECHNICAL  
SPECIFICATIONS

Technical Specifications...

<b>Dimensions:</b>		<b>Line Voltage:</b>	<b>Line Frequency:</b>
Height	15.30 cm (6.00 in.)	100-240 VAC	50 - 60 Hz
Width	31.10 cm (12.30 in.)	<b>Output Type</b>	<b>Duty Cycle</b>
Depth	41.30 cm (16.30 in.)		
Weight	< 8.8 kg (< 19 lbs)	RF ISOLATED <input checked="" type="checkbox"/> F	10 S/30 S
<b>Monopolar</b>	<b>Output Power</b>	<b>Output Frequency</b>	<b>Rep Rate</b>
Cut	300 W @ 300 Ω	490 kHz ± 5 kHz	N/A
Cut II	300 W @ 300 Ω	490 kHz ± 5 kHz	N/A
Blend (MAX)	200 W @ 300 Ω	490 kHz ± 5 kHz	30 kHz ± 5 kHz
Coagulation	120 W @ 500 Ω	490 kHz ± 5 kHz	30 kHz ± 5 kHz
Fulguration	80 W @ 500 Ω	490 kHz ± 5 kHz	25 kHz ± 5 kHz
Bipolar	80 W @ 150 Ω	490 kHz ± 5 kHz	30 kHz ± 5 kHz
<b>Warranty:</b> limited two (2) year on unit only			



Aaron® 3250™ Unit



BV-IDS-CS Stand  
(Optional)

Smoke Evacuation System  
(Optional)



Bovie Medical Corporation • 5115 Ulmerton Road • Clearwater, FL 33760 • www.boviemed.com • sales@boviemed.com  
U.S. Phone 1-800-537-2790 • Fax 1-800-323-1640 • Int'l Phone +1-727-384-2323 • Fax +1-727-347-9144

Bovie® and Aaron® are registered trademarks of Bovie Medical Corporation.